



ANTI-STATIC SPACER

CLEANING & INSTRUCTIONS FOR USE



Medical
Developments
International



Medical Developments International Limited
ABN 14 106 340 667

4 Caribbean Drive, Scoresby Victoria 3179, Australia
Telephone +61 3 9547 1888 Fax +61 3 9547 0262

Email mdi@medicaldev.com
Web www.medicaldev.com

8556MDI LI-BATAS-V1

INSTRUCTIONS FOR USE

<p>The Breath-A-Tech® Anti-static Spacer is a hand-held, portable valved holding chamber (VHC), that is placed between a metered dose inhaler (MDI, or "puffer") and the patient's mouth, to function as a reservoir from which an aerosol medication is dispensed in order to minimise delivery of large aerosolised particles. Intended for adults and children. The device is single patient use only. Can be used in most environments, including patient's homes, hospitals, emergency services and clinics.</p>					
<p>Remove the cap from your inhaler and shake well.</p>	<p>Push the inhaler into the base of the Breath-A-Tech® Anti-static Spacer</p>	<p>Put the mouthpiece between the lips and make an airtight seal, depress your inhaler, you will see a fine spray inside the Breath-A-Tech® Anti-static Spacer. Take a slow, deep breath through the mouth and hold for 5-10 seconds or continue to breathe in and out through the Breath-A-Tech® Anti-static Spacer for 4-6 breaths.</p>	<p>Wait at least one minute before repeating step 3. Only use one dose from the inhaler at a time.</p>	<p>If a face mask is needed, attach to the Breath-A-Tech® Anti-static Spacer. Firmly attach mask to the mouthpiece of the spacer and place mask over mouth and nose. Ensure the connection is secure at all times. Check for effective seal.</p>	<p>PLEASE NOTE Always follow the instructions provided with your asthma inhaler. The storage range is 5°C-40°C (41°F-104°F) at 15% to 95% relative humidity. ⚠ Before use check for foreign objects, broken or deteriorated parts. ⚠ Replace Breath-A-Tech® Anti-static Spacer immediately if any damage or missing parts are evident.</p>

CLEANING INSTRUCTIONS

<p>The Breath-A-Tech® Anti-static Spacer can be used straight out of the package, then should be cleaned at least once a month but can be washed more often if needed.</p>						
<p>The Breath-A-Tech® Anti-static Spacer can be used straight out of the package, then should be cleaned at least once a month but can be washed more often if needed.</p>	<p>Remove the base by placing your thumb in the hole then pull and twist at the same time. There is no need to dismantle the mouthpiece.</p>	<p>Wash all parts in warm water with a mild liquid detergent.</p>	<p>Rinse all parts in clean water.</p>	<p>Shake off excess water and allow to dry in a vertical position.</p>	<p>Firmly replace base once the Breath-A-Tech® Anti-static Spacer is dry.</p>	<p>DISHWASHING The Breath-A-Tech® Anti-static Spacer can be cleaned in a dishwasher up to 70°C. (up to 12 cycles). If dishwashing, follow step 1 and place the parts on top shelf of the dishwasher. Please place the spacer vertically in the dishwasher and ensure the outer side of the endcap faces upwards to prevent water pooling in device parts. If it is not clear to you on how to clean your Breath-A-Tech® Anti-static Spacer please ask your healthcare professional for help. REPLACE WITHIN 12 MONTHS</p>

SYMBOL		[LOT]	[REF]				
SYMBOL TITLE	Manufacturer	Batch code	Catalogue number	Consult instructions for use	Caution	Not made with natural rubber latex	Not made with BPA
EXPLANATORY TEXT	Indicates the medical device manufacturer.	Indicates the manufacturer's batch code so that the batch or lot can be identified.	Indicates the manufacturer's catalogue number so that the medical device can be identified.	Indicates the need for the user to use consult the instructions for use.	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	Not made with natural rubber latex	Not made with BPA